

(12) INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(19) World Intellectual Property
Organization
International Bureau



(43) International Publication Date
11 November 2004 (11.11.2004)

PCT

(10) International Publication Number
WO 2004/096430 A1

(51) International Patent Classification?: **B01J 19/02**

(21) International Application Number:
PCT/GB2004/001775

(22) International Filing Date: 26 April 2004 (26.04.2004)

(25) Filing Language: English

(26) Publication Language: English

(30) Priority Data:
0301260-6 29 April 2003 (29.04.2003) SE

(71) Applicant (for all designated States except MG, US): AS-
TRAZENECA AB [SE/SE]; Sodertalje, S-151 85 (SE).

(71) Applicant (for MG only): ASTRAZENECA UK LIM-
ITED [GB/GB]; 15 Stanhope Gate, London, Greater Lon-
don W1K 1LN (GB).

(72) Inventors; and

(75) Inventors/Applicants (for US only): BURNS, Stephen,
John [GB/GB]; AstraZeneca R & D Charnwood, Bakewell

Road, Loughborough, Leicestershire LE11 5RH (GB).
FLETCHER, Ian, Mark [GB/GB]; AstraZeneca R &
D Charnwood, Bakewell Road, Loughborough, Leices-
tershire LE11 5RH (GB). METCALF, Stephen, Paul
[GB/GB]; AstraZeneca R & D Charnwood, Bakewell
Road, Loughborough, Leicestershire LE11 5RH (GB).


(74) Agent: ASTRAZENECA; Global Intellectual Property,
S-151 85 Sodertalje (SE).

(81) Designated States (unless otherwise indicated, for every
kind of national protection available): AE, AG, AI, AM,
AT, AU, AZ, BA, BB, BG, BR, BW, BY, BZ, CA, CH, CN,
CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, EG, ES, FI,
GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE,
KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD,
MG, MK, MN, MW, MX, MZ, NA, NI, NO, NZ, OM, PG,
PH, PL, PT, RO, RU, SC, SD, SE, SG, SK, SL, SY, TJ, TM,
TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, YU, ZA, ZM,
ZW.

(84) Designated States (unless otherwise indicated, for every
kind of regional protection available): ARIPO (BW, GH,
GM, KE, LS, MW, MZ, NA, SD, SL, SZ, TZ, UG, ZM,
ZW), Eurasian (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM),
European (AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI,

[Continued on next page]

(54) Title: REACTOR COATED WITH A POLYMER



(57) Abstract: The invention relates to the use of a coated manufacturing vessel in pharmaceutical production. A polymer coating (ECTFE, PVDF, PFA) avoids the adhesion or sticking of compounds within the vessel.

WO 2004/096430 A1



FR, GB, GR, HU, IE, IT, LU, MC, NL, PL, PT, RO, SE, SI, SK, TR), OAPI (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG).

Published:

— with international search report

— before the expiration of the time limit for amending the claims and to be republished in the event of receipt of amendments

For two-letter codes and other abbreviations, refer to the "Guidance Notes on Codes and Abbreviations" appearing at the beginning of each regular issue of the PCT Gazette.

WO 2004/096430

PCT/GB2004/001775

1

REACTOR COATED WITH A POLYMER

5 The present invention relates to the use of a coated manufacturing vessel in pharmaceutical (MDI) production and its use with adhesive materials.

10 It is known that during the manufacture of ingredients of the pharmaceutical formulations can adhere to the surfaces of the manufacturing equipment. This is a potential problem when the formulation is dilute, the adhesive component is the active drug and therefore only a small amount of active drug substance is used. For example a drug such as formoterol (used as formoterol fumarate dihydrate) has a tendency to adhere to the surface of equipment used in its manufacture for pMDI's, and this results in significant losses of material (sometimes < 15%). The adhesion may be time dependent so that large scale manufacturing taking several days or where the batches are left stirring can result in
15 unpredictable deposition.

Overages of the active ingredient are allowed in limited amounts to compensate for this but as the adhesion is not reproducible, these overages may result in the variation of content of the active in the can and subsequent dosage of the drug from the pMDI.

20 In some circumstances the problem may be solved by manufacturing at sub-ambient temperature, and it is known in the art that optimising the manufacturing temperature may improve the product (Drug delivery to the Lungs, 1997, Production Scale Optimisation of the Manufacturing Process for HFA-134a metered Dose Inhaler Containing Both Salmeterol and Fluticasone propionate, SJ Duquemin et al) but this necessitates the use of
25 cooling equipment, the cost of which escalates as the scale of the manufacture increases. Also, using sub-ambient temperature manufacturing procedures has the possibility of affecting the behaviour of other ingredients, so sub-ambient procedures are best avoided.

30 The problem becomes even more significant if the ingredient is dilute as the losses constitute a greater proportion of the added amount than if there was a large amount added initially; also if the ingredient is expensive.

It has now been found that coating of the manufacturing equipment reduces losses of drug substance and negates the use of subambient temperatures for manufacture where one or

more of the ingredients of the mixture have a tendency to adhere to a conventional vessel surface.

In addition, cleaning of the equipment after use was found to be significantly easier than from an uncoated vessel where, due to the adhesion of ingredient to the uncoated surface, the cleaning was time-consuming. Efficient cleaning of pharmaceutical equipment is required under Good Manufacturing Practice.

In a first aspect the invention therefore provides use of a coated manufacturing vessel in the manufacture of a pharmaceutical formulation comprising a drug substance and pharmaceutically acceptable excipients.

Preferably the coated manufacturing vessel is used in the manufacture of a pharmaceutical formulation for a pMDI.

Preferably the coated manufacturing vessel is used in the manufacture of a pharmaceutical formulation for a pMDI, preferably wherein the active drug substance is formoterol fumarate dihydrate.

The vessels and processes of the invention are applicable to pressurised and non-pressurised manufacture e.g. aqueous suspensions for nebulisation, substances known to adhere to container surfaces esp steel e.g. proteins, antibodies.

As a further embodiment, the coating need not be limited to the vessel only, but stirrers, valves, tubing and any other equipment in the flow or in contact with the formulation mixture.

The vessel geometry was optimised for coating so that no voids were present where coating was difficult to apply and alterations were made to the standard design to prevent peeling of the coating from surfaces such as baffles and stirrer blades.

The coating can be any suitable alternative low surface energy coatings such as ECTFE (Ethylene Chloro-trifluoroethylene) PVDF (Polyvinylidene fluoride), and PFA (Perfluoroalkoxy) but the coating material is only limited by the suitability of the coating for manufacturing using pharmaceutical formulations (for example no leaching into, or contamination of, the formulation).

The manufacturing vessel is preferably coated with a PFA coating, preferably using a process described below. Preferably the vessel is a stainless steel vessel.

- 5 The invention is especially useful where the mixture is dilute and therefore the losses small but significant.

The following examples and data illustrate the invention.

10 Vessel Design and Coating Procedure

A suitable vessel is made from 316L Stainless steel, fabricated to British standard BS 5500 Category 2.

- 15 A suitable coating procedure is:

- Heat the vessel up to approx 400°C to burn off any manufacturing residues
- Apply a solvent based primer to the bare metal
- Heat the vessel to approx 400°C to cure the primer in place
- Coat the vessel with a thin layer of PFA polymer – the vessel is electrostatically charged, and the PFA applied as an inversely charged powder.
- 20 • Heat the vessel to approx 400°C to sinter the PFA onto the primer
- Repeat until a layer approx 0.5mm thick is build up.
- Test using a Holiday Detector – a commercially available tool for testing coating integrity. Measured conductivity between the vessel wall and the inside of the
- 25 vessel. Good coat = no electrical conduction

The invention applies to any design of vessel, one example of such a vessel, is represented by Figure 5

30 Process conditions for pMDI production using the coated vessel.

- A manufacturing process for Formoterol MDI HFA-134a:227 (75:25) blend was followed for the coated vessel. The coated vessel was used at a 20°C. No propellant headspace overage was required. Pressure filling was used, but the invention is not limited to
- 35 pressure filling.

The normal manufacturing procedure was followed except that the temperature of production was 20°C not the more usual 12°C. The lower temperature had been used previously as adhesion of components was reduced and the filling of and therefore dosing from, the canister was more reproducible.

Using the procedure above, a batch of pressurised inhalers was produced using the formulation given in Example 1 below:

Example 1

| Components | Quantity | |
|-------------------------------------------------------|----------|-----------------------|
| | % w/w | Theoretical per batch |
| Formoterol fumarate dihydrate micronised, conditioned | 0.0167 | 5.01 g |
| Polyethylene Glycol 1000 Ph. Eur. USP/NF | 0.10 | 30 g |
| Povidone (PVP) K-25 | 0.001 | 0.3 g |
| Propellant HFA-227 | 25 | 7.5 kg |
| Propellant HFA-134a | to 100 | to 30 kg |

All previous batches had been manufactured using an uncoated 30 Litre pressure vessel, some at 12°C, others at 20°C. Manufacturing time was normally 2 days, as was the case for this batch.

Product was filled into coated aerosol cans purged with propellant HFA-227, and crimped with a metering valve.

The results were compared with previous 30 Litre Giusti batches manufactured at ambient and sub-ambient temperatures.

Visual Observations

Observations were made through the sight glass on the vessel lid at various times throughout the batch. Deposition on the coated surfaces of the vessel side, stirrer and baffles was minimal. Three slight tidemarks were visible on the vessel side where the suspension had been left stirring overnight. No tidemarks were visible on the stirrer or baffles. The above observations were in contrast with batches made in a non-coated vessel, where tidemarks were clearly visible on all parts of the vessel.

Analysis of Results

The concentrations of Formoterol from samples taken at various points through the batch are given in Table 2.

Table 2: Formoterol Fumarate Dihydrate Concentration Through the Batch

| | | | | |
|--------------------------------|--------|--------|--------|--------|
| Amount of Suspension Used (kg) | 4.0 | 13.0 | 14.0 | 24.0 |
| Time after drug addition (hrs) | 18.5 | 20.5 | 43.0 | 45.5 |
| Mean* Concentration (%w/w) | 0.0163 | 0.0155 | 0.0154 | 0.0151 |

* based on a mean of three determinations

The above Table 2 shows that the concentration of active remains reasonably constant throughout the procedure.

If the mean concentration of Formoterol at any one point is divided by the theoretical concentration (based on weight added, in this case 0.0170% w/w) and multiplied by 100 then the % Formoterol concentration in the bulk mixture is given in Table 3.

Table 3: % Actual Formoterol Fumarate Dihydrate Concentration Through the Batch

| | | | | |
|--------------------------------|------|------|------|------|
| Amount of Suspension Used (kg) | 4.0 | 13.0 | 14.0 | 24.0 |
| Time after drug | 18.5 | 20.5 | 43.0 | 45.5 |

| addition (hrs) | | | | |
|-------------------------|-------|-------|-------|-------|
| Mean* Concentration (%) | 95.88 | 91.18 | 90.59 | 88.82 |

* of three cans per point.

These values are plotted in Figures 1 and 2.

5 The data shows a 4% loss of Formoterol immediately after component addition and a further 7% drop in the Formoterol concentration over 2 days.

However, significantly, negligible losses occurred when the mixture was held in the vessel for a significant time without processing. This is in contrast to previous batches and shows that the mixture can be left stirring overnight without loss in this type of vessel.

10 Also normally it is expected that the concentration of ingredients will increase as the level of mixture reduces due to propellant vapourising into the vessel headspace. This does not happen in the case of adhesive substances, because as the substance is deposited on the vessel walls the concentration decreases in the bulk. This process is unpredictable in the uncoated vessel but is more predicable in the coated vessel as can be seen in the graphical
15 representation.

This loss of ingredient is particularly significant where the amount is low ie dilute mixtures or suspensions as the loss on a large surface becomes more significant as the concentration of component decreases.

COMPARISON WITH OTHER BATCHES**Table 4: Losses of Active During the Manufacturing Process**

| Vessel Type | Temperature of Manufacture (°C) | Duration of Filling Period (days) | Loss of Formoterol at Start of Filling (%w/w) | Loss of Formoterol During Filling (%w/w) | Overall Loss of Formoterol at End of Batch (%w/w) |
|-------------|---------------------------------|-----------------------------------|-----------------------------------------------|------------------------------------------|---------------------------------------------------|
| Uncoated | 12 | 2 | 10.76 | 1.27 | 12.03 |
| Uncoated | 12 | 2 | 7.69 | 6.51 | 14.20 |
| Uncoated | 20 | 2 | 5.03 | 8.81 | 13.84 |
| Coated | 20 | 2 | 4.12 | 7.06 | 11.18 |

Table 4 shows that the coated vessel batch has the lowest overall loss of Formoterol at the end of the batch. The fall in concentration of Formoterol was found to be less than similar batches at 12 and 20C.

Interpretation of the Results

Using a manufacturing vessel all parts of which have been coated with a layer of material, which prevents any formulation component coming into the contact with any untreated surface, reduces the losses of that component onto the surfaces of the vessel.

It also results in a more predictable manufacturing process and is easier to clean than a corresponding uncoated vessel.

Figure 2 shows a Plot of % Formoterol Fumarate Dihydrate Concentration through the Batch against Amount of Suspension Used

Figure 3 shows a Plot of % Formoterol Fumarate Dihydrate Concentration compared to theoretical through the Batch against Time After Drug Addition

Figure 4: shows a smooth profile for Formoterol concentration through the coated vessel batch.

Figure 5: shows that the batch using the coated vessel has less loss at ambient compared
5 with batches in uncoated vessel at sub-ambient

Claims:

1. Use of a coated manufacturing vessel in the manufacture of a pharmaceutical formulation comprising a drug substance and pharmaceutically acceptable excipients.
2. Use according to claim 1 in the manufacture of a pharmaceutical formulation for a pMDI.
3. Use according to claim 1 or 2 in which the pharmaceutical formulation is a suspension for inhalation.
4. Use according to any one of claims 1 to 3 wherein the manufacture is carried out under pressurised conditions.
5. Use according to any one of claims 1 to 3 wherein the manufacture is carried out under non-pressurised conditions.
6. Use according to any one of claims 1 to 5 wherein the drug substance is formoterol fumarate dihydrate.
7. Use according to any one of claims 1 to 6 wherein the vessel is a stainless steel vessel.
8. Use according to any one of claims 1 to 7 wherein the coating is ECTFE, PVDF or PFA.
9. Use according to any one of claims 1 to 7 wherein the coating is PFA.
10. A manufacturing vessel characterised in that the vessel is coated to reduce adhesion of drug substance.
11. A stainless steel vessel according to claim 10 wherein the coating is ECTFE, PVDF or PFA.
12. A manufacturing vessel according to claim 10 wherein the coating is PFA.

13. A manufacturing vessel according to claim 11 or 12 which further includes coated piping, valves and accessories.

Figure 1




Figure 2




Figure 3




Figure 4





Figure 5



Coated regions of the vessel are shown in red (the thick black line in monotone- except the dip tube)

INTERNATIONAL SEARCH REPORT

International Application No.

PCT/GB2004/001775

A. CLASSIFICATION OF SUBJECT MATTER
IPC 7 B01J19/02

According to International Patent Classification (IPC) or to both national classification and IPC

B. FIELDS SEARCHED

Minimum documentation searched (classification system followed by classification symbols)

IPC 7 B01J

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched

Electronic data base consulted during the international search (name of data base and, where practical, search terms used)

EPO-Internal, WPI Data, PAJ

C. DOCUMENTS CONSIDERED TO BE RELEVANT

| Category * | Citation of document, with indication, where appropriate, of the relevant passages | Relevant to claim No. |
|------------|---------------------------------------------------------------------------------------------------------------------------------|-----------------------|
| X | EP 1 188 808 A (DAIKIN IND LTD). 20 March 2002 (2002-03-20) abstract page 2, paragraphs 2,5,13 page 3, paragraph 14 | 1-13 |
| X | DE 197 35 277 A (DEGUSSA) 18 February 1999 (1999-02-18) abstract page 2, line 1 - line 17 page 5; claims 1-4 | 1-13 |

☐ Further documents are listed in the continuation of box C.

☒ Parent family members are listed in annex.

* Special categories of cited documents:

- *A* document defining the general state of the art which is not considered to be of particular relevance
- *E* earlier document but published on or after the international filing date
- *L* document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)
- *O* document referring to an oral disclosure, use, exhibition or other means
- *P* document published prior to the international filing date but later than the priority date claimed

- *T* later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention
- *X* document of particular relevance; the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone
- *Y* document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art
- *Z* document member of the same patent family

Date of the actual completion of the international search

19 August 2004

Date of mailing of the international search report

01/09/2004

Name and mailing address of the ISA

European Patent Office, P.B. 5818 Patentlaan 2
NL - 2280 HV Rijswijk
Tel (+31-70) 340-2040, Tx. 31 651 epo nl,
Fax (+31-70) 340-3016

Authorized officer:

Thomasson, P

INTERNATIONAL SEARCH REPORT

International Application No
PCT/GB2004/001775

| Patent document cited in search report | | Publication date | | Patent family member(s) | Publication date |
|-------------------------------------------|---|---------------------|----|----------------------------|---------------------|
| EP 1188808 | A | 20-03-2002 | JP | 2000281939 A | 10-10-2000 |
| | | | EP | 1188808 A1 | 20-03-2002 |
| | | | US | 6734236 B1 | 11-05-2004 |
| | | | WO | 0058414 A1 | 05-10-2000 |
| DE 19735277 | A | 18-02-1999 | DE | 19735277 A1 | 18-02-1999 |